Fecal Immunochemical Test (FIT)

The FIT is a newer test that detects human hemoglobin in the stool. The test reacts to the protein component of the human hemoglobin and is more sensitive than guaiac-based tests that react to the heme component of the hemoglobin molecule for advanced adenomatous polyps and colon cancer.

Types of FIT

Manual Liquid-Based

Test results are read as either negative/positive Clinical Laboratory Improvement Amendments (CLIA)-waived Examples: Quidel QuickVue iFOB, Polymedco



Test results are read Hb in ng/ml or as negative/positive Not CLIA-waived

Examples: OC-Auto Micro 80, Polymedco

Manual Dry-Slide

Test results are read as either negative/positive

CLIA-waived


Examples: Beckman Coulter Hemoccult ICT, Immunostics Heme-Screen Specific





"The Perfect FIT"

- Detects globin rather than heme in hemoglobin, and therefore is highly specific for occult lower gastrointestinal (GI) bleeding as globin is largely degraded by upper GI enzymes
- No dietary or medication restrictions are required for stool collection
- Automated FIT has shown better performance characteristics than guaiac tests completed using three stool samples*
 - *Guaiac (Beckman Coulter Hemoccult II) sensitivity and specificity for detecting advanced adenomas (14%/92%), cancer (31%/92%), and advanced colorectal neoplasm (17%/93%)
 - *FIT (OC-Sensor) sensitivity and specificity for detecting advanced adenomas (34%/91%), cancer (85%/90%), and advanced colorectal neoplasm (44%/92%)1









More Information about FIT

Problems with Transporting FIT through the Mail

- Hemoglobin is not stable in stool samples
- Hemoglobin in stool samples stored in a liquid-based buffer degrades over time at a mean daily hemoglobin decrease of 29 ng Hg/ml²
- Prolonged delays in mailing or testing the sample result in more false-negative tests
- Mailed dry-slide stool samples have better stability of hemoglobin compared with liquid-based tests

Achieving the Best FIT Results

- The greater the number of stool samples provided, the higher the sensitivity for cancer or advanced colorectal neoplasia^{3,4,5} (Two consecutive stool samples are recommended)
- Patients should return or mail the stool sample on the day of collection
- The stool sample should be tested on the day of receipt

Important Points About Fecal Occult Blood Testing (FIT or Guaiac)

- FIT will not detect a polyp or tumor that is not bleeding
- Negative results may mean polyps are not bleeding
- Fecal occult blood tests should be completed annually
- Combining four randomized controlled trials, annual or biennial colorectal cancer (CRC) screening with fecal occult blood tests (FIT or guaiac), results indicated a significant reduction in colorectal cancer mortality of 16% by using the screening⁶
- Digital rectal exams should not be used to obtain stool for CRC screening
 - > Trauma of exam could cause bleeding
 - > For guaiac tests, sensitivities for detecting advanced neoplasia in 284 patients was 4.9% for digital rectal exam and 24% for a six-sample guaiac FOBT⁷
 - > For FIT, positive predictive value for CRC and large adenomatous polyps was 20% using a digital rectal exam and 27% in routine screening8
- A positive FIT indicates a need for a complete colonoscopy

References

¹Park D, Ryu S, Kim Y, et al., Comparison of guaiac-based and qualitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening. Am J Gastroenterol 2010;105:2017-2025.

²van Rossum LGM, van Rijn AF, van Oijen MGH, et al. False negative fecal occult blood tests due to delayed sample return in colorectal cancer screening. In. J. Cancer 2009;125:746-750.

³Yamamoto M, Nakama H. Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods. Hepato-Gastroenterology 2000;47:396-399.

⁴Levi Z, Rozen P, Hazazi R, et al. A quantitative immunochemical fecal occult blood test for colorectal neoplasia. Annals of Internal Medicine 2007;146:244-255. ⁵Park D, Ryu S, Kim YK, et al. Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening. Am J Gastroenterol 2010;105:2017-2025.

⁶Hewitson P, Glaziou P, Watson E, et al. Cochrane systematic review of colorectal cancer screening using the fecal occult blood test (Hemoccult): An update. Am J Gastroenterology 2008;103:1541-1549.

⁷Collins JF, Lieberman DA, Durbin TE, et al. Accuracy of screening for fecal occult blood on a single stool sample obtained by digital rectal examination: A comparison with recommended sampling practice. Ann Intern Med 2005;142:81-85.

8Nakama H, Zhang B. Does stool collection method affect outcomes in immunochemical fecal occult blood testing? Dis Colon Rectum 2001;44:871-875.







Collecting Stool Samples: A Resource for Patients

Collecting the Stool Sample

- Place a lined waste basket near your toilet
- Wash your hands
- Urinate before defecating. Then, flush the toilet to avoid getting urine on the stool sample
- Place rice paper on top of the water or seal the rice paper to the toilet to prevent contact with the water
- Rice paper may get wet
- Collect the specimen from the portion of the stool that is not in the water
- Use the probe or spatula to obtain the stool, but don't overfill the probe or spatula
- · Avoid overfilling the vial or card
- Discard spatula in waste basket
- Flush rice paper to discard

Do NOT Collect the Stool Sample When

- Menstruating
- · Bleeding hemorrhoids are present
- Blood is visible in the toilet
- · You have bleeding cuts on your hand
- · Toilet freshener is present in the toilet
- Rust or salt water are present in the toilet

When Collecting with Probe for a Vial

Collect from at least five different areas of the stool sample

When Collecting with Spatula for a Card

- · Collect from distinct areas as directed
- Mail or return sample to the office on the day of collection







